

Chlamydia Trachomatis HSP70 Recombinant Protein

Recombinant Chlamydia Trachomatis HSP70 (549-660)

CATALOG # AAA38097	HOST E Coli	APPLICATIONS ELISA
------------------------------	-----------------------	------------------------------

PRODUCT SYNONYMS

Chlamydia Trachomatis HSP70; Recombinant Chlamydia Trachomatis HSP70 (549-660); Chlamydia HSP70 (549-660a.a.); Chlamydia Trachomatis HSP70 (549-660 a.a.) Recombinant; Chlamydia HSP70; Chlamydia Trachomatis HSP70 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli
Specificity	Immunoreactive with sera of Chlamydia Trachomatis infected individuals.
Purity/Purification	Chlamydia HSP70 protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC. Chlamydia HSP70 protein was purified by proprietary chromatographic technique.
Form/Format	50mM Tris-HCl, pH 8.0, 8M Urea, 60mM NaCl and 50% glycerol.
Applicable Applications for Chlamydia Trachomatis HSP70 recombinant protein	ELISA
Preparation and Storage	Chlamydia HSP70 protein although stable at 4 degree C for 1 week, should be stored below -18 degree C. Please prevent freeze thaw cycles.

ADDITIONAL INFORMATION

Related Product Information for Chlamydia Trachomatis HSP70 recombinant protein	Description: The E Coli derived recombinant protein contains Chlamydia TrachomatisHSP70 protein epitopes, 549-660 amino acids. Introduction: Chlamydia is a common term for infection with any bacterium belonging to the phylum Chlamydiae. This term derives from the name of the bacterial genus Chlamydia in the family Chlamydiaceae, order Chlamydiales, class and phylum Chlamydiae. There are two genera in Chlamydiaceae: Chlamydia and Chlamydophila. The genus Chlamydia includes three species: C. trachomatis, C. muridarum, and C. suis.
Product Categories/Family for Chlamydia Trachomatis HSP70 recombinant protein	Viral Antigens; CHLAMYDIA TRACHOMATIS; Viral Antigens; Chlamydia

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/chlamydia-trachomatis-hsp70-recombinant-protein-e-coli-elisa/38097>