

Hemagglutinin-Influenza A Virus H3N2 Wisconsin 67/05 Recombinant Protein

H3N2 Influenza Virus Wisconsin/67/05 Recombinant H3N2 Wisconsin Recombinant

CATALOG #
AAA38380

PRODUCT SYNONYMS

Hemagglutinin-Influenza A Virus H3N2 Wisconsin 67/05; H3N2 Influenza Virus Wisconsin/67/05 Recombinant H3N2 Wisconsin Recombinant; Recombinant H3N2 Influenza Virus Wisconsin/67/05 H3N2 Wisconsin Recombinant; Hemagglutinin-Influenza A Virus H3N2 Wisconsin 67/05 recombinant protein

PRODUCT SPECIFICATIONS

Purity/Purification	Greater than 90.0% as determined by SDS-PAGE.
Form/Format	The Recombinant H3N2 A/Wisconsin/67/05 solution (715ug/ml) contains 10mM Sodium phosphate, pH 7.2 and 150mM Sodium Chloride. and 0.005% Tween-20
Source	Baculovirus Insect Cells
Physical Appearance	Sterile Filtered colorless solution
Preparation and Storage	H3N2 A/Wisconsin/67/05 Recombinant should be stored at 4 degree C. Do not freeze!

ADDITIONAL INFORMATION

Related Product Information for Hemagglutinin-Influenza A Virus H3N2 Wisconsin 67/05 recombinant protein	Description: Recombinant Full-Length H3N2 A/Wisconsin/67/05 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells and its Mw is 72,000 dalton. Introduction: H3N2 is a subtype of the influenza A virus. Its name derives from the forms of the two kinds of protein on the surface of its coat, hemagglutinin(H) and neuraminidase(N). H3N2 exchanges genes for internal proteins with other influenza subtypes. H3N2 has tended to dominate in prevalence over H1N1, H1N2, and influenza B. H3N2 strain descended from H2N2 by antigenic shift, in which genes from multiple subtypes re-assorted to form a new virus. Both the H2N2 and H3N2 strains contained genes from avian influenza viruses.
Product Categories/Family for Hemagglutinin-Influenza A Virus H3N2 Wisconsin 67/05 recombinant protein	Viral Antigens; INFLUENZA VIRUS; Viral Antigens; Influenza

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/hemagglutinin-influenza-a-virus-h3n2-wisconsin-67/05-recombinant-protein/38380>