

SARS Associated Spike Mosaic S(N) Recombinant Protein | SARS recombinant protein

Recombinant SARS Associated Spike Mosaic S(N)

CATALOG # AAA38534	HOST E.Coli	APPLICATIONS Western Blot, ELISA
------------------------------	-----------------------	--

PRODUCT SYNONYMS

SARS Associated Spike Mosaic S(N); Recombinant SARS Associated Spike Mosaic S(N); SARS Mosaic S(N); SARS-Associated Coronavirus Spike Mosaic S (N) Recombinant; SARS S(N); SARS recombinant protein

PRODUCT SPECIFICATIONS

Host	E.Coli
Specificity	SARS Mosaic protein is Immunoreactive with sera of SARS-infected individuals.
Purity/Purification	SARS Mosaic Protein is >95% pure as determined by 10% PAGE (coomassie staining). SARS Mosaic protein was purified by proprietary chromatographic technique.
Form/Format	(1mg/ml) 25mM Tris-HCl, 0.4% sarcosyl, 0.25% Triton 100 and 50% glycerol.
Applicable Applications for SARS recombinant protein	WB (Western Blot), ELISA
Physical Appearance	Sterile Filtered clear solution.
Preparation and Storage	SARS Mosaic protein is shipped at ambient temperature. Upon arrival, store at -20C.

ADDITIONAL INFORMATION

Related Product Information for SARS recombinant protein	Description: The E Coli derived 38 kDa mosaic protein contains the N-terminal section of the Spike protein 12-53, 90-115, 171-203 amino acids immunodominant regions. Introduction: SARS Coronavirus is an enveloped virus containing three outer structural proteins, namely the membrane (M), envelope (E), and spike (S) proteins. Spike (S)-glycoprotein of the virus interacts with a cellular receptor and mediates membrane fusion to allow viral entry into susceptible target cells. Accordingly, S-protein plays an important role in virus infection cycle and is the primary target of neutralizing antibodies.
Product Categories/Family for SARS recombinant protein	Viral Antigens; SARS PROTEINS; Viral Antigens; SARS

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/sars-associated-spike-mosaic-s-n-recombinant-protein-e-coli-wb-western-blot-elisa/38534>