

Porcine P-selectin, sP-Selectin ELISA Kit

Porcine soluble P-selectin, sP-Selectin ELISA Kit

CATALOG #
AAA41266

REACTIVITY
Porcine


PRODUCT SYNONYMS

P-selectin; sP-Selectin; Porcine soluble P-selectin; sP-Selectin ELISA Kit; sP-Selectin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.2ng/ml - 60ng/ml
Sensitivity	0.12ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for P-selectin, sP-Selectin elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Porcine Soluble E-selectin (also known as SE-S) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Porcine SE-S antibody. SE-S present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Porcine SE-S Antibody is added and binds to SE-S in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated SE-S antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Porcine SE-S. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/p-selectin-sp-selectin/41266>