

Human Toll-like receptor 5, TLR5 ELISA Kit | TLR5 elisa kit

Human Toll-like receptor 5, TLR5 ELISA Kit

CATALOG #
AAA41634

REACTIVITY
Human


PRODUCT SYNONYMS

Toll-like receptor 5; TLR5; Human Toll-like receptor 5; TLR5 ELISA Kit; TLR5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	871
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.1ng/ml - 40ng/ml
Sensitivity	0.05ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TLR5 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of human Toll-like Receptor 5 (also known as TLR5) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with human TLR5 antibody. TLR5 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated human TLR5 Antibody is added and binds to TLR5 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated TLR5 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of human TLR5. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	224967052
NCBI GeneID	289337
NCBI Accession #	NP_001139300.1
NCBI GenBank Nucleotide #	NM_001145828.1
NCBI Official Full Name	toll-like receptor 5
NCBI Official Synonym Full Names	toll-like receptor 5
NCBI Official Symbol	Tlr5
NCBI Protein Information	toll-like receptor 5

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/toll-like-receptor-5-tlr5/41634>