

Human matrix metalloproteinase 5, MMP-5 ELISA Kit | MMP-5 elisa kit

Human matrix metalloproteinase 5, MMP-5 ELISA Kit

CATALOG #
AAA42437

REACTIVITY
Human


PRODUCT SYNONYMS

matrix metalloproteinase 5; MMP-5; Human matrix metalloproteinase 5; MMP-5 ELISA Kit; MMP-5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.1ng/ml - 30ng/ml
Sensitivity	0.06ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MMP-5 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of human Matrix Metalloproteinase 5 (also known as MMP-5) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with human MMP-5 antibody. MMP-5 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated human MMP-5 Antibody is added and binds to MMP-5 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated MMP-5 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of human MMP-5. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/matrix-metalloproteinase-5-mmp-5/42437>