

Mouse Renalase (RNLS) ELISA Kit | RNLS elisa kit

Mouse Renalase (RNLS) ELISA Kit

CATALOG #
AAA43495

REACTIVITY
Mouse


PRODUCT SYNONYMS

Renalase (RNLS); Mouse Renalase (RNLS) ELISA Kit; RNLS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	51
Samples	Serum, plasma, cell culture supernates, ascites, tissue homogenates or other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	0.5-100U/L
Sensitivity	0.249U/L
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for RNLS elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Mouse Renalase (also known as RNLS) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Mouse RNLS antibody. RNLS present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Mouse RNLS Antibody is added and binds to RNLS in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated RNLS antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Mouse RNLS. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	363902006
Molecular Weight	5,372 Da
NCBI Official Full Name	renalase, partial
UniProt Protein Name	Renalase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/renalase-rnls/43495>