

Human mid-regional pro-atrial natriuretic peptide (MR-proANP) ELISA Kit | MR-proANP elisa kit

Human mid-regional pro-atrial natriuretic peptide (MR-proANP) ELISA Kit

CATALOG # AAA43856	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

mid-regional pro-atrial natriuretic peptide (MR-proANP); Human mid-regional pro-atrial natriuretic peptide (MR-proANP) ELISA Kit; MR-proANP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	5-2000pmol/L
Sensitivity	2.35pmol/L
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MR-proANP elisa kit	Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human MR-proANP antibody. MR-proANP present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human MR-proANP Antibody is added and binds to MR-proANP in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated MR-proANP antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human MR-proANP. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm. Intended Uses: This sandwich kit is for the accurate quantitative detection of Human mid-regional pro-atrial natriuretic peptide (also known as MR-proANP) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.
---	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mid-regional-pro-atrial-natriuretic-peptide-mr-proanp/43856>