

Human Mitochondrial import inner membrane translocase subunit Tim8 A (TIMM8A) ELISA Kit | TIMM8A elisa kit

Human Mitochondrial import inner membrane translocase subunit Tim8 A (TIMM8A) ELISA Kit

CATALOG # AAA44342	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

Mitochondrial import inner membrane translocase subunit Tim8 A (TIMM8A); Human Mitochondrial import inner membrane translocase subunit Tim8 A (TIMM8A) ELISA Kit; TIMM8A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	20ng/L - 3500ng/L
Sensitivity	9.56ng/L
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TIMM8A elisa kit	Intended Uses: This sandwich kit is for the accurate quantitative detection of human Mitochondrial Import Inner Membrane Translocase Subunit Tim8 A (also known as TIMM8A) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates. Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with human TIMM8A antibody. TIMM8A present in the sample is added and binds to antibodies coated on the wells. And then biotinylated human TIMM8A Antibody is added and binds to TIMM8A in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated TIMM8A antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of human TIMM8A. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mitochondrial-import-inner-membrane-translocase-subunit-tim8-a-timm8a/44342>