

Cy3 Conjugated Goat Anti-Rabbit IgG Secondary Antibody

Cy3 Conjugated Goat Anti-Rabbit IgG secondary antibody

CATALOG # AAA45614	HOST Goat	APPLICATIONS Immunocytochemistry, Immunohistochemistry, Immunohistochemistry
------------------------------	---------------------	--


PRODUCT SYNONYMS

Cy3 Conjugated Goat Anti-Rabbit IgG; Polyclonal Secondary Antibody; Cy3 Conjugated Goat Anti-Rabbit IgG secondary antibody

PRODUCT SPECIFICATIONS


Host	Goat
Clonality	Polyclonal
Specificity	This Cy3 conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with mouse, bovine IgG.
Form/Format	Concentrated, Liquid
Concentration	2mg/ml (varies by lot)
Applicable Applications for Cy3 Conjugated Goat Anti-Rabbit IgG secondary antibody	ICC (Immunocytochemistry), IHC (Immunohistochemistry), IHC (Immunohistochemistry)
Preparation and Storage	At 4 degree C for one year.

IMAGES



IHC (Immunohistochemistry)

Collagen IV was detected in paraffin-embedded sections of rat kidney tissues using rabbit anti-Collagen IV Antigen Affinity purified polyclonal antibody at 1 ug/mL. Cy3 Conjugated Goat Anti-Rabbit IgG secondary antibody was used to detect the primary antibody at 20ug/mL.



IHC (Immunohistochemistry)

MCM2 was detected in paraffin-embedded sections of human intestinal cancer tissues using rabbit anti-MCM2 Antigen Affinity purified polyclonal antibody at 1ug/mL. Cy3 Conjugated Goat Anti-Rabbit IgG secondary antibody was used to detect the primary antibody at 20ug/mL.

ADDITIONAL INFORMATION

Related Product Information for Cy3 Conjugated Goat Anti-Rabbit IgG secondary antibody

Description: This antibody is purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody for rabbit IgG. The antibody preparation is solid phase adsorbed with human serum proteins and bovine serum proteins to ensure minimal cross reactivity in tissue or cell preparations.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

