

Humanized anti-Human CLDN6 Monoclonal Antibody | anti-CLDN6 antibody

Anti-CLDN6 (IMAB027) mAb

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA47573	Humanized	Human	Flow Cytometry, Functional Assay, ELISA


PRODUCT SYNONYMS

CLDN6; Monoclonal Antibody; Anti-CLDN6 (IMAB027) mAb; Claudin 6; Claudin-6; Skullin; anti-CLDN6 antibody

PRODUCT SPECIFICATIONS


Host	Humanized
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG1
Purity/Purification	Purified from cell culture supernatant by affinity chromatography
Form/Format	Powder. Lyophilized from sterile PBS, pH7.4. Normally 5%-8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Applicable Applications for anti-CLDN6 antibody	FCM/FACS (Flow Cytometry), ELISA
Preparation and Storage	Store at -20 degree C to -80 degree C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80 degree C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

IMAGES




FCM/FACS (Flow Cytometry)

Figure 1. Flow cytometry analysis with 1 ug/mL Anti-CLDN6 (IMAB027) mAb (AAA47573) on Expi293 cells transfected with Human CLDN6 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



FCM/FACS (Flow Cytometry)

Figure 1. Flow cytometry analysis with 1 ug/mL Anti-CLDN6 (IMAB027) mAb (AAA47573) on Expi293 cells transfected with Human CLDN6 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



ELISA
Figure 1. Flow cytometry analysis with 1 ug/mL Anti-CLDN6 (IMAB027) mAb (AAA47573) on Expi293 cells transfected with Human CLDN6 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

ADDITIONAL INFORMATION

Product Categories/Family for anti-CLDN6 antibody	Biosimilar Benchmark Antibodies
---	---------------------------------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/monoclonal/cldn6-antibody-flow-cytometry-elisa/47573>