

Mouse PROGESTERONE Monoclonal Antibody

MOUSE ANTI PROGESTERONE

CATALOG # AAA49204	HOST Mouse	APPLICATIONS immunoassay, Radioimmunoassay, ELISA
------------------------------	----------------------	---

PRODUCT SYNONYMS

PROGESTERONE; Monoclonal Antibody; MOUSE ANTI PROGESTERONE; anti-PROGESTERONE antibody

PRODUCT SPECIFICATIONS

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone Number	BGN/6-5E-10B
Specificity	Mouse anti Progesterone antibody, clone BGN/6-5E-10B has the following specificities: Progesterone-11a-HMS-BSA 100% 5-beta-Pregnane-3,20-dione 48% 5-alpha-Pregnane-3,20-dione 26.4% 17-alpha-Hydroxyprogesterone 2.5% 20-alpha-Hydroxyprogesterone 0.04%
Form/Format	Purified Purified IgG - Liquid
Concentration	IgG concentration 1.0 mg/ml (varies by lot)
Applicable Applications for anti-PROGESTERONE antibody	RIA (Radioimmunoassay), ELISA
Perservative Stabilisers	0.09% Sodium Azide
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Immunogen	Progesterone-11a-hemisuccinate conjugated to bovine serum albumin (Pr11a-HMS-BSA)
Buffer Solution	Phosphate buffered saline
Target Species	Broad
Preparation and Storage	Store at 4 degree C or at -20 degree C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. Shelf Life: 18 months from date of despatch.

ADDITIONAL INFORMATION

Related Product Information for anti-PROGESTERONE antibody	Mouse anti Progesterone antibody, clone BGN/6-5E-10B has the following specificities: Progesterone-11a-HMS-BSA100%5-beta-Pregnane-3,20-dione48%5-alpha-Pregnane-3,20-dione26.4%17-alpha-Hydroxyprogesterone2.5%20-alpha-Hydroxyprogesterone0.04%
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

