

POLY CASPASE Assay Kit

SR FLICA POLY CASPASE KIT

CATALOG #
AAA50646

APPLICATIONS
Immunofluorescence, Flow Cytometry, Functional Assay


PRODUCT SYNONYMS

POLY CASPASE; SR FLICA POLY CASPASE KIT; POLY CASPASE assay kit

PRODUCT SPECIFICATIONS

Form/Format	Kit
Reagents In The Kits	1 vial of SR-VAD-FMK FLICA Reagent - lyophilized 10X Wash Buffer, 15 mL Fixative, 6 mL Hoechst Stain, 1 mL
Applicable Applications for POLY CASPASE assay kit	IF (Immunofluorescence), FCM/FACS (Flow Cytometry)
Preparation and Storage	Store the unopened kit (and each unopened component) at 4 degree C until the expiration date. Protect the FLICA reagent from light at all times. Once reconstituted, the 150X FLICA stock should be stored at -20 degree C protected from light. Shelf Life: The 150X FLICA reconstituted stock is stable for up to 6 months and may be thawed twice during that time.

IMAGES



Application Data

Staining of staurosporine-induced Jurkat cells using the SR FLICATM poly caspase kit

ADDITIONAL INFORMATION

Related Product Information for POLY CASPASE assay kit	Description: The SR FLICA Poly Caspase kit uses a quick and easy method to analyze active caspases in apoptotic cells. This FLICA reagent binds the following activated caspases: Caspase 1, Caspase 3, Caspase 4, Caspase 5, Caspase 6, Caspase 7, Caspase 8 and Caspase 9. It can also be used as a generic probe for the detection of most caspases. Test Principle: Caspase FLICA kits measure apoptosis by detecting active caspases in whole, living cells. These kits do not work by using antibodies or as an ELISA. Instead, their methodology is based on a unique cell-permeable and non-cytotoxic reagent called the Fluorochrome Inhibitor of Caspases (FLICA). The FLICA reagent contains a caspase inhibitor sequence linked to a red (Sulforhodamine, SR) fluorescent probe. The Caspase FLICA Kits are suitable for cells in suspension, adherent cells, thin tissue sections (but not fixed or paraffin-embedded cells) from many species including mammalian, insect and yeast. Different cell types, e.g. HeLa, primary neurons, macrophages and lymphocytes have also been successfully studied with these kits. This kit can be used with a fluorescence microscope or a fluorescence plate reader using black microtitre plates.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/poly-caspase/50646>