

Human Coagulation Factor X (F10) ELISA Kit | F10 elisa kit

Human Coagulation Factor X (F10) ELISA Kit

CATALOG #
AAA52800

REACTIVITY
Human


PRODUCT SYNONYMS

Coagulation Factor X (F10); Human Coagulation Factor X (F10) ELISA Kit; F10 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	488
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for F10 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human F10 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

182390

NCBI GeneID	2159
Molecular Weight	54,732 Da
NCBI Official Full Name	coagulation factor X
NCBI Official Synonym Full Names	coagulation factor X
NCBI Official Symbol	F10
NCBI Official Synonym Symbols	FX; FXA
NCBI Protein Information	coagulation factor X; factor Xa; prothrombinase; Stuart-Prower factor
UniProt Protein Name	Coagulation factor X
UniProt Gene Name	F10
UniProt Entry Name	FA10_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/coagulation-factor-x-f10/52800>