

Porcine Toll Like Receptor 9 (TLR-9) ELISA Kit | TLR-9 elisa kit

Porcine Toll Like Receptor 9 (TLR-9) ELISA Kit

CATALOG #
AAA52955

REACTIVITY
Porcine


PRODUCT SYNONYMS

Toll Like Receptor 9 (TLR-9); Porcine Toll Like Receptor 9 (TLR-9) ELISA Kit; TLR-9 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Sequence Length	1032
Samples	Porcine serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant TLR-9 concentration
Assay Type	Sandwich
Detection Range	20 ng/ml-0.312 ng/ml
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TLR-9 elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Porcine TLR-9 monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	157057166
NCBI GeneID	81897
NCBI Accession #	NP_112455.2
NCBI GenBank Nucleotide #	NM_031178.2

Molecular Weight	116,412 Da
NCBI Official Full Name	toll-like receptor 9
NCBI Official Synonym Full Names	toll-like receptor 9
NCBI Official Symbol	Tlr9
NCBI Protein Information	toll-like receptor 9
UniProt Protein Name	Toll-like receptor 9
UniProt Gene Name	Tlr9
UniProt Entry Name	TLR9_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/toll-like-receptor-9-tlr-9/52955>