

Mouse PSA ELISA Kit | PSA elisa kit

Mouse PSA ELISA Kit

CATALOG #
AAA53133

REACTIVITY
Mouse


PRODUCT SYNONYMS

PSA; Mouse PSA ELISA Kit; PSA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	261
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSA elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Mouse PSA monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4502173
NCBI GeneID	354
NCBI Accession #	NP_001639.1

NCBI GenBank Nucleotide #	NM_001648.2
Molecular Weight	25,116 Da
NCBI Official Full Name	prostate-specific antigen isoform 1 preproprotein
NCBI Official Synonym Full Names	kallikrein-related peptidase 3
NCBI Official Symbol	KLK3
NCBI Official Synonym Symbols	APS; PSA; hK3; KLK2A1
NCBI Protein Information	prostate-specific antigen; seminin; P-30 antigen; kallikrein-3; semenogelase; gamma-seminoprotein; prostate specific antigen
UniProt Protein Name	Prostate-specific antigen
UniProt Gene Name	KLK3
UniProt Synonym Gene Names	APS; PSA; Semin
UniProt Entry Name	KLK3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/psa/53133>