

Human Protamine 1 (PRM1) ELISA Kit | PRM1 elisa kit

Human Protamine 1 (PRM1) ELISA Kit

CATALOG #
AAA53186

REACTIVITY
Human


PRODUCT SYNONYMS

Protamine 1 (PRM1); Human Protamine 1 (PRM1) ELISA Kit; PRM1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	51
Assay Type	Sandwich
Detection Range	1000 pg/mL-15.6 pg/mL
Sensitivity	Up to 5 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PRM1 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human ANP monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

119605563

NCBI GeneID	5619
Molecular Weight	6,823 Da
NCBI Official Full Name	protamine 1
NCBI Official Synonym Full Names	protamine 1
NCBI Official Symbol	PRM1
NCBI Official Synonym Symbols	P1; CT94.1
NCBI Protein Information	sperm protamine P1; cysteine-rich protamine; testis specific protamine 1; cancer/testis antigen family 94, member 1
UniProt Protein Name	Sperm protamine P1
UniProt Gene Name	PRM1
UniProt Entry Name	HSP1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/protamine-1-prm1/53186>