

Porcine Lactoperoxidase (LPO) ELISA Kit | LPO elisa kit

Porcine Lactoperoxidase (LPO) ELISA Kit

CATALOG #
AAA53527

REACTIVITY
Porcine


PRODUCT SYNONYMS

Lactoperoxidase (LPO); Porcine Lactoperoxidase (LPO) ELISA Kit; LPO elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Specificity	No cross-reaction with other factors.
Sequence Length	324
Samples	Serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	10 ng/ml-0.156 ng/ml
Sensitivity	Up to 0.05 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LPO elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Porcine LPO monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	187224
NCBI GeneID	4025

Molecular Weight	70,984 Da
NCBI Official Full Name	lactoperoxidase, partial
NCBI Official Synonym Full Names	lactoperoxidase
NCBI Official Symbol	LPO
NCBI Official Synonym Symbols	SPO
NCBI Protein Information	lactoperoxidase; salivary peroxidase
UniProt Protein Name	Lactoperoxidase
UniProt Gene Name	LPO
UniProt Synonym Gene Names	SAPX; LPO; SPO
UniProt Entry Name	PERL_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/lactoperoxidase-lpo/53527>