

Human Adenosine Diphosphate (ADP) ELISA Kit | ADP elisa kit

Human Adenosine Diphosphate (ADP) ELISA Kit

CATALOG #
AAA53745

REACTIVITY
Human


PRODUCT SYNONYMS

Adenosine Diphosphate (ADP); Human Adenosine Diphosphate (ADP) ELISA Kit; ADP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	628
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADP elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human ADP monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

665401814

NCBI GeneID	37073
NCBI Accession #	NP_001286563.1
NCBI GenBank Nucleotide #	NM_001299634.1
NCBI Official Full Name	adipose, isoform B
NCBI Official Synonym Full Names	adipose
NCBI Official Symbol	adp
NCBI Official Synonym Symbols	anon-WO0196371.1; anon-WO0196371.27; anon-WO0196371.3; CG5124; DmelCG5124
NCBI Protein Information	CG5124 gene product from transcript CG5124-RA; CG5124-PA; CG5124-PB; adipose; adp-PA; adp-PB

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/adenosine-diphosphate-adp/53745>