

Human hyaluronidase (HAase) ELISA Kit | HAase elisa kit

Human hyaluronidase (HAase) ELISA Kit

CATALOG #
AAA53778

REACTIVITY
Human


PRODUCT SYNONYMS

hyaluronidase (HAase); Human hyaluronidase (HAase) ELISA Kit; HAase elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	463
Samples	Human serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant HAase concentration
Assay Type	Sandwich
Detection Range	200 ng/ml-3.12 ng/ml
Sensitivity	0.6 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HAase elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human HAase monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	2935328
NCBI GeneID	8372

Molecular Weight	15,865 Da
NCBI Official Full Name	hyaluronidase
NCBI Official Synonym Full Names	hyaluronoglucosaminidase 3
NCBI Official Symbol	HYAL3
NCBI Official Synonym Symbols	LUCA3; HYAL-3; LUCA-3
NCBI Protein Information	hyaluronidase-3; lung carcinoma protein 3
UniProt Protein Name	Hyaluronidase-3
UniProt Gene Name	HYAL3
UniProt Synonym Gene Names	LUCA3; Hyal-3; LuCa-3
UniProt Entry Name	HYAL3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/hyaluronidase-haase/53778>