

Porcine MPO-ANCA IgG ELISA Kit | MPO-ANCA IgG elisa kit

Porcine MPO-ANCA IgG ELISA Kit

CATALOG # AAA53830	REACTIVITY Porcine
------------------------------	------------------------------


PRODUCT SYNONYMS

MPO-ANCA IgG; Porcine MPO-ANCA IgG ELISA Kit; MPO-ANCA IgG elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Sequence Length	745
Assay Type	Sandwich
Samples	Porcine serum, plasma, or cell culture supernatant
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MPO-ANCA IgG elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antigen Sandwich ELISA technique. The pre-coated antigen is a related-Porcine MPO-ANCA IgG antigen, while the detection antigen is another antigen. Samples are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4557759
NCBI GeneID	4353

NCBI Accession #	NP_000241.1
NCBI GenBank Nucleotide #	NM_000250.1
Molecular Weight	87,248 Da
NCBI Official Full Name	myeloperoxidase
NCBI Official Synonym Full Names	myeloperoxidase
NCBI Official Symbol	MPO
NCBI Protein Information	myeloperoxidase
UniProt Protein Name	Myeloperoxidase
UniProt Gene Name	MPO
UniProt Synonym Gene Names	MPO
UniProt Entry Name	PERM_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/mpo-anca-igg/53830>