

Rat Deoxyribonuclease I (DNASE1) ELISA Kit | DNASE1 elisa kit

Rat Deoxyribonuclease I (DNASE1) ELISA Kit

CATALOG # AAA54085	REACTIVITY Rat
------------------------------	--------------------------


PRODUCT SYNONYMS

Deoxyribonuclease I (DNASE1); Rat Deoxyribonuclease I (DNASE1) ELISA Kit; DNASE1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	No cross-reaction with other factors.
Sequence Length	282
Assay Type	Sandwich
Detection Range	100 ng/mL-1.56 ng/mL
Sensitivity	Up to 0.5 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DNASE1 elisa kit	Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Rat DNASE1 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	1197173
Molecular Weight	31,434 Da

NCBI Official Full Name	deoxyribonuclease I
UniProt Protein Name	Deoxyribonuclease-1
UniProt Gene Name	DNASE1
UniProt Synonym Gene Names	DNL1; DRNI; DNase I
UniProt Entry Name	DNAS1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/deoxyribonuclease-i-dnase1/54085>