

Human Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit | CKMB elisa kit

Human Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit

CATALOG #
AAA54162

REACTIVITY
Human


PRODUCT SYNONYMS

Creatine Kinase MB Isoenzyme (CKMB); Human Creatine Kinase MB Isoenzyme (CKMB) ELISA Kit; CKMB elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	381
Samples	Human serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	Up to 0.5 ng/ml
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CKMB elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human CKMB monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	21536288
NCBI GeneID	1158

NCBI Accession #	NP_001815.2
NCBI GenBank Nucleotide #	NM_001824.4
Molecular Weight	43,101 Da
NCBI Official Full Name	creatine kinase M-type
NCBI Official Synonym Full Names	creatine kinase, muscle
NCBI Official Symbol	CKM
NCBI Official Synonym Symbols	CKMM; M-CK
NCBI Protein Information	creatine kinase M-type; creatine kinase M chain
UniProt Protein Name	Creatine kinase M-type
UniProt Gene Name	CKM
UniProt Synonym Gene Names	CKMM
UniProt Entry Name	KCRM_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/creatine-kinase-mb-isoenzyme-ckmb/54162>