

Porcine Dopamine Receptor D1 (DRD1) ELISA Kit | DRD1 elisa kit

Porcine Dopamine Receptor D1 (DRD1) ELISA Kit

CATALOG # AAA54553	REACTIVITY Porcine
------------------------------	------------------------------


PRODUCT SYNONYMS

Dopamine Receptor D1 (DRD1); Porcine Dopamine Receptor D1 (DRD1) ELISA Kit; DRD1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Specificity	No cross-reaction with other factors.
Sequence Length	446
Samples	Serum, plasma or cell culture supernatant and organizations in the natural and recombinant DRD1 concentration.
Assay Type	Sandwich
Detection Range	10 ng/ml-0.156 ng/ml
Sensitivity	0.05 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DRD1 elisa kit	Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Porcine DRD1 monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	20269979
NCBI GeneID	1812

Molecular Weight	49,293 Da
NCBI Official Full Name	dopamine receptor D1
NCBI Official Synonym Full Names	dopamine receptor D1
NCBI Official Symbol	DRD1
NCBI Official Synonym Symbols	DADR; DRD1A
NCBI Protein Information	D(1A) dopamine receptor; D(1A) dopamine receptor; dopamine D1 receptor
UniProt Protein Name	D(1A) dopamine receptor
UniProt Gene Name	DRD1
UniProt Entry Name	DRD1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/dopamine-receptor-d1-drd1/54553>