

Human filamentous hemagglutinin (FHA) ELISA Kit | FHA elisa kit

Human filamentous hemagglutinin (FHA) ELISA Kit

CATALOG #
AAA54689

REACTIVITY
Human


PRODUCT SYNONYMS

filamentous hemagglutinin (FHA); Human filamentous hemagglutinin (FHA) ELISA Kit; FHA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	29202
Samples	Human serum, plasma or cell culture supernatant and organizations in the natural and recombinant FHA concentration.
Assay Type	Sandwich
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	0.5 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FHA elisa kit

Principle of the Assay||This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human FHA monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	334143384
NCBI GeneID	10719085

NCBI Accession #	YP_004536540.1
NCBI GenBank Nucleotide #	NC_015581.1
Molecular Weight	2,894,531 Da
NCBI Official Full Name	filamentous hemagglutinin family outer membrane protein
NCBI Official Symbol	Thicy_0285
NCBI Protein Information	filamentous hemagglutinin
UniProt Protein Name	Filamentous hemagglutinin family outer membrane protein
UniProt Entry Name	F6DA95_THICA

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/filamentous-hemagglutinin-fha/54689>