

Human immunosuppressive acidic protein (IAP) ELISA Kit | IAP elisa kit

Human immunosuppressive acidic protein (IAP) ELISA Kit

CATALOG #
AAA54816

REACTIVITY
Human


PRODUCT SYNONYMS

immunosuppressive acidic protein (IAP); Human immunosuppressive acidic protein (IAP) ELISA Kit; IAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	224
Samples	Human serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	Up to 0.5 ng/ml
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IAP elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human IAP monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	33622285
NCBI GeneID	1725101
NCBI Accession #	NP_891932.1

NCBI GenBank Nucleotide #	NC_005068.1
Molecular Weight	25,829 Da
NCBI Official Full Name	iap
NCBI Official Symbol	ClgVgp085
NCBI Protein Information	orf85; iap
UniProt Protein Name	Iap
UniProt Entry Name	Q7T5K8_GVCL

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/immunosuppressive-acidic-protein-iap/54816>