

Chicken Bone alkaline phosphatase (BALP) ELISA Kit | BALP elisa kit

Chicken Bone alkaline phosphatase (BALP) ELISA Kit

CATALOG #
AAA54878

REACTIVITY
Chicken


PRODUCT SYNONYMS

Bone alkaline phosphatase (BALP); Chicken Bone alkaline phosphatase (BALP) ELISA Kit; BALP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Chicken
Specificity	No cross-reaction with other factors.
Samples	Chicken serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant BALP concentration
Assay Type	Sandwich
Detection Range	200 U/L-3.12 U/L
Sensitivity	the minimum detectable Chicken BALP up to 0.6 U/L.
Intra-assay Precision	? 8%
Inter-assay Precision	? 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BALP elisa kit

Principle of the assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Chicken BALP monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.