

Chicken Advanced oxidation protein products (AOPP) ELISA Kit | AOPP elisa kit

Chicken Advanced oxidation protein products (AOPP) ELISA Kit

CATALOG #
AAA54981

REACTIVITY
Chicken


PRODUCT SYNONYMS

Advanced oxidation protein products (AOPP); Chicken Advanced oxidation protein products (AOPP) ELISA Kit; AOPP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Chicken
Sequence Length	555
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AOPP elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Chicken AOPP monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	31542077
NCBI GeneID	103213

NCBI Accession #	NP_598761.2
NCBI GenBank Nucleotide #	NM_134000.3
NCBI Official Full Name	adapter protein CIKS
NCBI Official Synonym Full Names	TRAF3 interacting protein 2
NCBI Official Symbol	Traf3ip2
NCBI Official Synonym Symbols	Act1; CIKS; AI429613
NCBI Protein Information	adapter protein CIKS

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/chicken-elisa-kits/advanced-oxidation-protein-products-aopp/54981>