

Human Follistatin (FS) ELISA Kit | FS elisa kit

Human Follistatin (FS) ELISA Kit

CATALOG #
AAA55019

REACTIVITY
Human


PRODUCT SYNONYMS

Follistatin (FS); Human Follistatin (FS) ELISA Kit; FS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	317
Samples	Human serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant FS concentration
Assay Type	Sandwich
Detection Range	200 ng/ml-3.12 ng/ml
Sensitivity	0.6 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FS elisa kit

Principle of the assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human FS monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	182721
NCBI GeneID	10468

Molecular Weight	34,803 Da
NCBI Official Full Name	follicle-stimulating hormone receptor
NCBI Official Synonym Full Names	follicle-stimulating hormone receptor
NCBI Official Symbol	FST
NCBI Official Synonym Symbols	FS
NCBI Protein Information	follicle-stimulating hormone receptor; activin-binding protein; follicle-stimulating hormone receptor isoform FST317
UniProt Protein Name	Follicle-stimulating hormone receptor
UniProt Gene Name	FST
UniProt Synonym Gene Names	FS
UniProt Entry Name	FST_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/follicle-stimulating-hormone-receptor-fs/55019>