

Porcine Luteinizing Hormone (LH) ELISA Kit | LH elisa kit

Porcine Luteinizing Hormone (LH) ELISA Kit

CATALOG #
AAA55045

REACTIVITY
Porcine


PRODUCT SYNONYMS

Luteinizing Hormone (LH); Porcine Luteinizing Hormone (LH) ELISA Kit; LH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Sequence Length	5
Samples	Porcine serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant LH concentration
Assay Type	Sandwich
Detection Range	50 mIU/ml-0.78 mIU/ml
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LH elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Porcine LH monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4262812
NCBI GeneID	3972
Molecular Weight	15,345 Da
NCBI Official Full Name	luteinizing hormone beta subunit (LHbeta), partial
NCBI Official Synonym Full Names	luteinizing hormone beta polypeptide

NCBI Official Symbol	LHB
NCBI Official Synonym Symbols	CGB4; hLHB; LSH-B
NCBI Protein Information	lutropin subunit beta; LSH-beta; lutropin beta chain; luteinizing hormone beta subunit; interstitial cell stimulating hormone, beta chain
UniProt Protein Name	Lutropin subunit beta
UniProt Gene Name	LHB
UniProt Synonym Gene Names	Lutropin beta chain; LH-B; LSH-B; LSH-beta
UniProt Entry Name	LSHB_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/luteinizing-hormone-lh/55045>