

Porcine Relaxin (RLN) ELISA Kit | RLN elisa kit

Porcine Relaxin (RLN) ELISA Kit

CATALOG #
AAA55108

REACTIVITY
Porcine


PRODUCT SYNONYMS

Relaxin (RLN); Porcine Relaxin (RLN) ELISA Kit; RLN elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Specificity	No cross-reaction with other factors.
Sequence Length	185
Samples	Porcine serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant RLN concentration
Assay Type	Sandwich
Detection Range	2000 pg/ml-31.2 pg/ml
Sensitivity	12 pg/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for RLN elisa kit

Principle of the assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Porcine RLN monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	490056
NCBI GeneID	6013

Molecular Weight	13,011 Da
NCBI Official Full Name	relaxin
NCBI Official Synonym Full Names	relaxin 1
NCBI Official Symbol	RLN1
NCBI Official Synonym Symbols	H1; H1RLX; RLXH1; bA12D24.3.1; bA12D24.3.2
NCBI Protein Information	prorelaxin H1; preprorelaxin H1
UniProt Protein Name	Prorelaxin H1
UniProt Gene Name	RLN1
UniProt Entry Name	REL1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/relaxin-rln/55108>