

Porcine Thrombomodulin (TM) ELISA Kit | TM elisa kit

Porcine Thrombomodulin (TM) ELISA Kit

CATALOG #
AAA55135

REACTIVITY
Porcine


PRODUCT SYNONYMS

Thrombomodulin (TM); Porcine Thrombomodulin (TM) ELISA Kit; TM elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Sequence Length	575
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TM elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Porcine TM monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507483
NCBI GeneID	7056
NCBI Accession #	NP_000352.1

NCBI GenBank Nucleotide #	NM_000361.2
Molecular Weight	60,329 Da
NCBI Official Full Name	thrombomodulin
NCBI Official Synonym Full Names	thrombomodulin
NCBI Official Symbol	THBD
NCBI Official Synonym Symbols	TM; THRM; AHUS6; BDCA3; CD141; THPH12
NCBI Protein Information	thrombomodulin; fetomodulin; CD141 antigen
UniProt Protein Name	Thrombomodulin
UniProt Gene Name	THBD
UniProt Synonym Gene Names	THRM; TM
UniProt Entry Name	TRBM_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/thrombomodulin-tm/55135>