

Mouse Bovine Serum Albumin (BSA) ELISA Kit | BSA elisa kit

Mouse Bovine Serum Albumin (BSA) ELISA Kit

CATALOG # AAA55200	REACTIVITY Mouse
------------------------------	----------------------------


PRODUCT SYNONYMS

Bovine Serum Albumin (BSA); Mouse Bovine Serum Albumin (BSA) ELISA Kit; BSA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	607
Samples	Mouse serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant BSA concentration
Assay Type	Sandwich
Detection Range	200 ug/ml-3.12 ug/ml
Sensitivity	0.6 ug/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BSA elisa kit	Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Mouse BSA monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	3336842
NCBI GeneID	280717

UniProt Accession #	P02769; O02787; P04277; Q3SZR2; A5PJX3
Molecular Weight	69,293 Da
NCBI Official Full Name	bovine serum albumin
NCBI Official Synonym Full Names	albumin
NCBI Official Symbol	ALB
NCBI Protein Information	serum albumin; BSA
UniProt Protein Name	Serum albumin
UniProt Gene Name	ALB
UniProt Entry Name	ALBU_BOVIN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/bovine-serum-albumin-bsa/55200>