

Porcine Sorbitol Dehydrogenase (SDH) ELISA Kit | SDH elisa kit

Porcine Sorbitol Dehydrogenase (SDH) ELISA Kit

CATALOG # AAA55203	REACTIVITY Porcine
------------------------------	------------------------------


PRODUCT SYNONYMS

Sorbitol Dehydrogenase (SDH); Porcine Sorbitol Dehydrogenase (SDH) ELISA Kit; SDH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Specificity	No cross-reaction with other factors.
Sequence Length	357
Samples	Porcine serum, plasma or cell culture supernatant
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	0.5 ng/mL
Intra-assay Precision	? 8%
Inter-assay Precision	? 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SDH elisa kit	Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Porcine SDH monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	156627571
NCBI GeneID	6652
NCBI Accession #	NP_003095.2

NCBI GenBank Nucleotide #	NM_003104.5
Molecular Weight	38,325 Da
NCBI Official Full Name	sorbitol dehydrogenase
NCBI Official Synonym Full Names	sorbitol dehydrogenase
NCBI Official Symbol	SORD
NCBI Official Synonym Symbols	SORD1; HEL-S-95n
NCBI Protein Information	sorbitol dehydrogenase; L-iditol 2-dehydrogenase; epididymis secretory sperm binding protein Li 95n
UniProt Protein Name	Sorbitol dehydrogenase
UniProt Gene Name	SORD
UniProt Entry Name	DHSO_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/sorbitol-dehydrogenase-sdh/55203>