

Human Collagen Type IV (COL4) ELISA Kit | COL4 elisa kit

Human Collagen Type IV (COL4) Elisa Kit

CATALOG #
AAA55273

REACTIVITY
Human


PRODUCT SYNONYMS

Collagen Type IV (COL4); Human Collagen Type IV (COL4) Elisa Kit; COL4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	18
Samples	Human serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant COL4 concentration
Assay Type	Sandwich
Detection Range	10 ng/ml-0.156 ng/ml
Sensitivity	0.05 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for COL4 elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human COL4 monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4261609
NCBI GeneID	1287

Molecular Weight	161,632 Da
NCBI Official Full Name	collagen type IV alpha 5 chain, partial
NCBI Official Synonym Full Names	collagen, type IV, alpha 5
NCBI Official Symbol	COL4A5
NCBI Official Synonym Symbols	ATS; ASLN; CA54
NCBI Protein Information	collagen alpha-5(IV) chain; collagen alpha-5(IV) chain; dA24A23.1; dA149D17.3; collagen IV, alpha-5 polypeptide; collagen of basement membrane, alpha-5 chain
UniProt Protein Name	Collagen alpha-5(IV) chain
UniProt Gene Name	COL4A5
UniProt Entry Name	CO4A5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/collagen-type-iv-col4/55273>