

Human campylobacter jejuni (CJ) antibody (IgM) ELISA Kit | CJ elisa kit

Human campylobacter jejuni (CJ) antibody (IgM) ELISA Kit

CATALOG #
AAA55283

REACTIVITY
Human


PRODUCT SYNONYMS

campylobacter jejuni (CJ) antibody (IgM); Human campylobacter jejuni (CJ) antibody (IgM) ELISA Kit; CJ elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	194
Samples	Human serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant CJ antibody (IgM) concentration
Assay Type	Sandwich
Detection Range	100 U/L-1.56U/L
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CJ elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human CJ antibody (IgM) monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	166007338
NCBI Official Full Name	Chain A, Crystal Structure Of Apo-Pgld From Campylobacter Jejuni

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.