

## Dog Ferritin heavy chain ELISA Kit | FTH1 elisa kit

Dog Ferritin heavy chain ELISA Kit

CATALOG #  
**AAA57334**

REACTIVITY  
**Dog**


### PRODUCT SYNONYMS

Ferritin heavy chain; Dog Ferritin heavy chain ELISA Kit; Ferritin H subunit; FTH1; FTH; 1.16.3.1; FTH1 elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                     |
|-------------------------|---------------------------------------------------------------------|
| Reactivity              | Dog                                                                 |
| Sequence Length         | 183                                                                 |
| Detection Range         | 0.625-40 ng/mL                                                      |
| Preparation and Storage | For long term storage, please store the entire kit at -20 degree C. |

### IMAGES



Standard Curve (Sample)

### NCBI AND UNIPROT INFORMATION

|                           |                      |
|---------------------------|----------------------|
| NCBI GI #                 | 50978756             |
| NCBI GeneID               | 403631               |
| NCBI Accession #          | NP_001003080.1       |
| NCBI GenBank Nucleotide # | NM_001003080.1       |
| UniProt Accession #       | Q95MP7; Q53VC0       |
| Molecular Weight          | 21,308 Da            |
| NCBI Official Full Name   | ferritin heavy chain |
| NCBI Official Symbol      | FTH1                 |
| NCBI Protein Information  | ferritin heavy chain |
| UniProt Protein Name      | Ferritin heavy chain |
| UniProt Gene Name         | <b>FTH1</b>          |

|                            |                         |
|----------------------------|-------------------------|
| UniProt Synonym Gene Names | FTH; Ferritin H subunit |
| UniProt Entry Name         | FRIH_CANLF              |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/fth1-dog-elisa-kits/ferritin-heavy-chain/57334>