

Clostridium difficile Toxin A Monoclonal Antibody

MAb to C. difficile Toxin A

CATALOG # AAA57665	HOST Host: Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from Balb/c mice.Source: Ascites	APPLICATIONS ELISA
------------------------------	---	------------------------------

PRODUCT SYNONYMS

Clostridium difficile Toxin A; Monoclonal Antibody; MAb to C. difficile Toxin A; Monoclonal Antibody to Clostridium difficile Toxin A; anti-Clostridium difficile Toxin A antibody

PRODUCT SPECIFICATIONS

Host	Host: Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from Balb/c mice. Source: Ascites
Clonality	Monoclonal
Isotype	IgG2a
Specificity	Clostridium difficile Toxin A C. difficile Toxin A. Does not cross react with C. difficile Toxin B.
Form/Format	Ascites, Liquid
Concentration	6.44mg/ml (estimated antibody concentration determined by densitometry scan) (varies by lot)
Applicable Applications for anti-Clostridium difficile Toxin A antibody	ELISA
Immunogen	C. difficile Toxin A
Affinity Constant	Not determined
Buffer	Not applicable
Preservative	0.1% Sodium Azide
Dry Ice Note	The product may be shipped with dry ice, and an additional fee may be added to your shipping cost.
Important Note	Centrifuge before opening to ensure complete recovery of vial contents.
Warnings	This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.
Preparation and Storage	Upon receipt, store at -40 degree C. Avoid multiple freeze/thaw cycles.

ADDITIONAL INFORMATION

Product Categories/Family for anti-Clostridium difficile Toxin A antibody	Monoclonal Antibodies to Infectious Agents and Toxins
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.