

Soluble Tumor Necrosis Factor Receptor 2 (STNF-R2) Recombinant Protein | STNF-R2 recombinant protein

Soluble Tumor Necrosis Factor Receptor 2 (STNF-R2)

CATALOG #

AAA58768


PRODUCT SYNONYMS

Soluble Tumor Necrosis Factor Receptor 2 (STNF-R2); sTNFR II; Tumor receptor beta receptor; Tumor necrosis factor binding protein 2; Tumor necrosis factor receptor superfamily member 1B; Tumor necrosis factor receptor type II; p80 TNF alpha receptor; Etanercept; CD120b; p75; TBPII; TNFRII; TNF-R75; TNFBR; TNFR2; TNFR80.; STNF-R2 recombinant protein

PRODUCT SPECIFICATIONS

Purity/Purification	>95%
Form/Format	20mM Tris-HCl pH8.0, 5mM EDTA, 50% Glycerol
Concentration	0.30mg/ml (varies by lot)
Preparation and Storage	Store -20°C. Vials should be centrifuge briefly before opening to ensure contents are not lost when removing cap)

IMAGES



SDS-PAGE

ADDITIONAL INFORMATION

Related Product Information for STNF-R2 recombinant protein	Two types of soluble TNF receptors (sTNFR-I and sTNFR-II) have been identified, which act to neutralize the biological activities of TNFalpha and TNF beta. Levels of these soluble receptors appear to arise as a result of shedding of the extracellular domains of the membrane bound receptors. Elevated levels of soluble TNF receptors have been found in the amniotic fluid of pregnant women. Recombinant sTNFR-I comprises a 161 amino acid fragment (41-201) corresponding to the mature sTNFR-I protein and is expressed in E Coli with an amino-terminal hexahistidine tag. This product is not intended for therapeutic use.
Product Categories/Family for STNF-R2 recombinant protein	Recombinant Proteins; Soluble Tumor Necrosis Factor Receptor 2
References	Wallach, D. et al., (1991) Agents Actions Suppl. 35: 51- 57

NCBI AND UNIPROT INFORMATION

Molecular Weight	24.45kDa
------------------	----------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.