

Human Progesterone Receptor ELISA Kit | PGR elisa kit

Progesterone Receptor Human ELISA Kit

CATALOG # AAA58895	REACTIVITY Human	APPLICATIONS ELISA
------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Progesterone Receptor; Progesterone Receptor Human ELISA Kit; Human Progesterone Receptor; PGR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Rat Serum, Plasma, Urine, Cell Culture Supernatant Or Tissue Samples
Assay Type	Sandwich
Detection Range	31.2 pg/ml-2000 pg/ml
Applicable Applications for PGR elisa kit	ELISA
Preparation and Storage	Store at 2-8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGR elisa kit

Principle of the Assay: This kit was based on standard sandwich enzyme-linked immune-sorbent assay technology. The purified anti-KIM-1 antibody was pre-coated onto 96-well plates. And the HRP conjugated anti-KIM-1 antibody was used as detection antibodies. The standards, test samples and HRP conjugated detection antibody were added to the wells subsequently, mixed and incubated, then, unbound conjugates were washed away with wash buffer. TMB substrates (A & B) were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the KIM-1 amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of KIM-1 can be calculated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	329112574
NCBI GeneID	100848466
NCBI Accession #	NP_001192285.1
NCBI GenBank Nucleotide #	NM_001205356.1
NCBI Official Full Name	progesterone receptor
NCBI Official Symbol	PGR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/progesterone-receptor/58895>