

Nuclear and Cytoplasmic Protein Extraction Kit

Nuclear and Cytoplasmic Protein Extraction Kit

CATALOG #
AAA59049

APPLICATIONS
Western Blot

PRODUCT SYNONYMS

Nuclear and Cytoplasmic Protein Extraction Kit; Nuclear and Cytoplasmic Protein Extraction Kit extraction kit

PRODUCT SPECIFICATIONS

Applicable Applications for Nuclear and Cytoplasmic Protein Extraction Kit extraction kit

WB (Western Blot)

ADDITIONAL INFORMATION

Related Product Information for Nuclear and Cytoplasmic Protein Extraction Kit extraction kit

1. Nuclear and Cytoplasmic Protein Extraction Kit is designed to quickly prepare highly enriched fractions of cytoplasmic and nuclear proteins from eukaryotic samples, such as cultured cells and fresh tissues.
2. This kit provides a simple and convenient method for extracting nuclear and cytoplasmic protein. The user can complete the separation of nuclear and cytoplasmic protein within about 90 min. Extracted proteins can be directly apply for downstream reactions, such as Western, EMSA, foot printing, reporter gene detection, and measurement of enzyme activity.

Product Categories/Family for Nuclear and Cytoplasmic Protein Extraction Kit extraction kit

[Lysis & Protein Extraction](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/extraction-kits/nuclear-and-cytoplasmic-protein-extraction-kit/59049>