

Rabbit OLIG2 Monoclonal Antibody | anti-OLIG2 antibody

OLIG2 (RBT-OLIG2), RMaB, 0.1mL Concentrate

CATALOG # AAA59472	HOST Rabbit	REACTIVITY Paraffin, Frozen
------------------------------	-----------------------	---------------------------------------

PRODUCT SYNONYMS

OLIG2; Monoclonal Antibody; OLIG2 (RBT-OLIG2); RMaB; 0.1mL Concentrate; (RBT-OLIG2); anti-OLIG2; anti OLIG2; Oligodendrocytes transcription factor 2; anti-Oligodendrocytes transcription factor 2; anti Oligodendrocytes transcription factor 2; OLIG-2 oligodendrocyte transcription factor 2; basic helix-loop-helix (bHLH) transcription factor; BHLHB1; motor neuron and oligodendrocyte differentiation modulator; Serine (Ser) Ser30; Ser10; Ser13; Ser14 phosphoprotein; RACK17; anti-OLIG-2 oligodendrocyte transcription factor 2; anti OLIG-2 oligodendrocyte transcription factor 2; anti-basic helix-loop-helix (bHLH) transcription factor; anti basic helix-loop-helix (bHLH) transcription factor; anti-BHLHB1; anti BHLHB1; anti-motor neuron and oligodendrocyte differentiation modulator; anti motor neuron and oligodendrocyte differentiation modulator; anti-Serine (Ser) Ser30; anti Serine (Ser) Ser30; anti-RACK17; anti RACK17; H10; H 10; anti-OLIG2 antibody


PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Paraffin, Frozen
Clonality	Monoclonal
Isotype	IgG
Clone Number	RBT-OLIG2
Control	Brain Astrocytoma

IMAGES

Antibody Type	Dilution	Volume/Conc
Prediluted	Ready-To-Use	3.0 ml
Prediluted	Ready-To-Use	7.0 ml
Prediluted	Ready-To-Use	15.0 ml
Concentrated	1:50-1:200	0.1 ml
Concentrated	1:50-1:200	0.5 ml
Concentrated	1:50-1:200	1.0 ml
Control Slides		5

Dilution Info



IHC (Immunohistochemistry)

IHC of OLIG2 Antibody on a FFPE Neuroblastoma Tissue

ADDITIONAL INFORMATION

Product Categories/Family for anti-OLIG2 antibody	Nuclear; Cytoplasmic
---	----------------------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.