

## Rapid-Lenti HIV-1 p24 ELISA Kit | HIV-1 elisa kit

Rapid-Lenti HIV-1 p24 ELISA Kit (Single Wash)

CATALOG #  
**AAA60734**

### PRODUCT SYNONYMS


Rapid-Lenti HIV-1 p24; Rapid-Lenti HIV-1 p24 ELISA Kit (Single Wash); HIV-1 elisa kit

### PRODUCT SPECIFICATIONS

Preparation and Storage

Store at +4 degree C

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for HIV-1 elisa kit

Principle of the Assay: Microtitration wells coated with murine anti-HIV-1 p24 capture antibody are exposed to test specimens, which may contain HIV-1 p24 reactive determinants. A mixture of biotinylated anti-HIV-1 p24 antibodies and streptavidin HRP conjugate are then added forming a complex (HIV-1 p24/antibodies/conjugate) that binds to the wells. After an incubation period, unbound components in the test sample are washed away. Specifically bound enzyme conjugate is detected by reaction with the Substrate Solution, tetramethylbenzidine (TMB). The assay is measured spectrophotometrically to indicate the level of HIV-1 p24 reactive determinants present in a sample. The assay can be used to monitor HIV-1 in cell culture or to determine the viral titer of lentiviral samples and can be completed in as little as 90 minutes.

#### Product Categories/Family for HIV-1 elisa kit

**HIV RESEARCH PRODUCTS**

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/elisa-kits/rapid-lenti-hiv-1-p24/60734>