

Human Adenylyl Cyclase Associated Protein 2 (CAP2) ELISA Kit | CAP2 elisa kit

Human Adenylyl Cyclase Associated Protein 2 (CAP2) ELISA Kit

CATALOG #
AAA64299

REACTIVITY
Human


PRODUCT SYNONYMS

Adenylyl Cyclase Associated Protein 2 (CAP2); Human Adenylyl Cyclase Associated Protein 2 (CAP2) ELISA Kit; CAP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of CAP2. No significant cross-reactivity or interference between CAP2 and analogues was observed.
Sequence Length	477
Samples	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.156-10ng/mL
Sensitivity	0.056ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level CAP2 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level CAP2 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CAP2 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of CAP2 in human serum, plasma, tissue homogenates, cell lysates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to CAP2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to CAP2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain CAP2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of CAP2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	5453593
NCBI GeneID	10486
NCBI Accession #	NP_006357.1
NCBI GenBank Nucleotide #	NM_006366.2
UniProt Accession #	P40123; Q6IAY2; B2R5Y3; B7Z1C4; B7Z214
Molecular Weight	24,018 Da
NCBI Official Full Name	adenylyl cyclase-associated protein 2
NCBI Official Synonym Full Names	CAP, adenylate cyclase-associated protein, 2 (yeast)
NCBI Official Symbol	CAP2
NCBI Protein Information	adenylyl cyclase-associated protein 2
UniProt Protein Name	Adenylyl cyclase-associated protein 2
UniProt Gene Name	CAP2
UniProt Synonym Gene Names	CAP 2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.