

Human Proteasome Subunit Alpha Type 1 (PSMa1) ELISA Kit | PSMa1 elisa kit

Human Proteasome Subunit Alpha Type 1 (PSMa1) ELISA Kit

CATALOG #
AAA64306

REACTIVITY
Human


PRODUCT SYNONYMS

Proteasome Subunit Alpha Type 1 (PSMa1); Human Proteasome Subunit Alpha Type 1 (PSMa1) ELISA Kit; HC2; NU; PROS30; Proteasome(Prosome; Macropain)Subunit; Alpha Type; 1; Macropain subunit C2; Multicatalytic endopeptidase complex subunit C2; Proteasome nu chain; PSMa1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of PSMa1. No significant cross-reactivity or interference between PSMa1 and analogues was observed.
Sequence Length	263
Samples	Tissue homogenates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.156-10ng/mL
Sensitivity	0.067ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level PSMa1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level PSMa1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSMA1 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of PSMA1 in human tissue homogenates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PSMA1. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to PSMA1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PSMA1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of PSMA1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4506179
NCBI GeneID	5682
NCBI Accession #	NP_002777.1
NCBI GenBank Nucleotide #	NM_002786.3
UniProt Accession #	P25786; Q53YE8; Q9BRV9; A8K400
Molecular Weight	30,239 Da
NCBI Official Full Name	proteasome subunit alpha type-1 isoform 2
NCBI Official Synonym Full Names	proteasome subunit alpha 1
NCBI Official Symbol	PSMA1
NCBI Official Synonym Symbols	NU; HC2; PROS30; HEL-S-275
NCBI Protein Information	proteasome subunit alpha type-1
UniProt Protein Name	Proteasome subunit alpha type-1
UniProt Gene Name	PSMA1
UniProt Synonym Gene Names	HC2; NU; PROS30; PSC2; PROS-30

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/proteasome-subunit-alpha-type-1-psma1/64306>