

Human Bcl2 Associated Athanogene 3 (BAG3) ELISA Kit | BAG3 elisa kit

Human Bcl2 Associated Athanogene 3 (BAG3) ELISA Kit

CATALOG #
AAA64768

REACTIVITY
Human


PRODUCT SYNONYMS

Bcl2 Associated Athanogene 3 (BAG3); Human Bcl2 Associated Athanogene 3 (BAG3) ELISA Kit; BIS; CAIR1; BAG Family Molecular Chaperone Regulator 3; Bcl-2-binding protein Bis; Docking protein CAIR-1; BAG3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of BAG3. No significant cross-reactivity or interference between BAG3 and analogues was observed.
Sequence Length	575
Samples	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.156-10ng/mL
Sensitivity	0.069ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level BAG3 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level BAG3 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BAG3 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of BAG3 in human serum, plasma, tissue homogenates, cell lysates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to BAG3. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to BAG3. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain BAG3, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of BAG3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	14043024
NCBI GeneID	9531
NCBI Accession #	NP_004272.2
NCBI GenBank Nucleotide #	NM_004281.3
UniProt Accession #	O95817; Q3B763; Q9NT20; Q9P120; A8K5L8
Molecular Weight	61,595 Da
NCBI Official Full Name	BAG family molecular chaperone regulator 3
NCBI Official Synonym Full Names	BCL2 associated athanogene 3
NCBI Official Symbol	BAG3
NCBI Official Synonym Symbols	BIS; MFM6; BAG-3; CAIR-1
NCBI Protein Information	BAG family molecular chaperone regulator 3
UniProt Protein Name	BAG family molecular chaperone regulator 3
UniProt Gene Name	BAG3
UniProt Synonym Gene Names	BIS; BAG-3

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/bcl2-associated-athanogene-3-bag3/64768>