

Human Stress Induced Phosphoprotein 1 (STIP1) ELISA Kit | STIP1 elisa kit

Human Stress Induced Phosphoprotein 1 (STIP1) ELISA Kit

CATALOG #
AAA65062

REACTIVITY
Human


PRODUCT SYNONYMS

Stress Induced Phosphoprotein 1 (STIP1); Human Stress Induced Phosphoprotein 1 (STIP1) ELISA Kit; HOP; P60; IEF-SSP-3521; STI1L; Hsp70/Hsp90-Organizing Protein; Renal carcinoma antigen NY-REN-11; Transformation-sensitive protein IEF SSP 3521; STIP1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of STIP1. No significant cross-reactivity or interference between STIP1 and analogues was observed.
Sequence Length	590
Samples	Serum, plasma, tissue homogenates and other biological fluids.
Assay Type	Sandwich
Detection Range	15.625-1000pg/mL
Sensitivity	6.5pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level STIP1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level STIP1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for STIP1 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of STIP1 in human serum, plasma, tissue homogenates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to STIP1. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to STIP1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain STIP1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of STIP1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	544063423
NCBI GeneID	10963
NCBI Accession #	NP_001269581.1
NCBI GenBank Nucleotide #	NM_001282652.1
UniProt Accession #	P31948; Q3ZCU9; Q5TZU9; B4DM70; F5H0T1; G3XAD8
Molecular Weight	59,778 Da
NCBI Official Full Name	stress-induced-phosphoprotein 1 isoform a
NCBI Official Synonym Full Names	stress induced phosphoprotein 1
NCBI Official Symbol	STIP1
NCBI Official Synonym Symbols	HOP; P60; STI1; STI1L; HEL-S-94n; IEF-SSP-3521
NCBI Protein Information	stress-induced-phosphoprotein 1
UniProt Protein Name	Stress-induced-phosphoprotein 1
UniProt Gene Name	STIP1
UniProt Synonym Gene Names	STI1; Hop

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.