

Human Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit | TICAM1 elisa kit

Human Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit

CATALOG #
AAA65290

REACTIVITY
Human


PRODUCT SYNONYMS

Toll Like Receptor Adaptor Molecule 1 (TICAM1); Human Toll Like Receptor Adaptor Molecule 1 (TICAM1) ELISA Kit; NAK; T2K; TICAM1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of TICAM1. No significant cross-reactivity or interference between TICAM1 and analogues was observed.
Sequence Length	712
Samples	Tissue homogenates, cell lysates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.312-20ng/mL
Sensitivity	0.124ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level TICAM1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level TICAM1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TICAM1 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of TICAM1 in human tissue homogenates, cell lysates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to TICAM1. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to TICAM1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain TICAM1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of TICAM1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	41281981
NCBI GeneID	148022
NCBI Accession #	NP_891549.1
NCBI GenBank Nucleotide #	NM_182919.3
UniProt Accession #	Q8IUC6; O75532; Q86XP8; Q96GA0; B3Y691
Molecular Weight	76,422 Da
NCBI Official Full Name	TIR domain-containing adapter molecule 1
NCBI Official Synonym Full Names	toll like receptor adaptor molecule 1
NCBI Official Symbol	TICAM1
NCBI Official Synonym Symbols	TRIF; IIAE6; MyD88-3; PRVTIRB; TICAM-1
NCBI Protein Information	TIR domain-containing adapter molecule 1
UniProt Protein Name	TIR domain-containing adapter molecule 1
UniProt Gene Name	TICAM1
UniProt Synonym Gene Names	PRVTIRB; TRIF; TICAM-1; MyD88-3; TIR domain-containing adapter protein inducing IFN-beta

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/toll-like-receptor-adaptor-molecule-1-ticam1/65290>