

Human N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit | NT-ProBNP elisa kit

Human N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit

CATALOG # AAA65475	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP); Human N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit; PAF-Acether; AGEPC; Acetyl-Glycerol-Ether-Phosphorylcholine; NT-ProBNP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Assay Type	Sandwich
Detection Range	39.062-2500pg/mL
Sensitivity	15pg/mL
Intra Assay	CV<10%
Inter Assay	CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/n-terminal-pro-brain-natriuretic-peptide-nt-probnp/65475>