

Human Paraoxonase 2 (PON2) ELISA Kit | PON2 elisa kit

Human Paraoxonase 2 (PON2) ELISA Kit

CATALOG #
AAA66720

REACTIVITY
Human


PRODUCT SYNONYMS

Paraoxonase 2 (PON2); Human Paraoxonase 2 (PON2) ELISA Kit; Serum Paraoxonase/Arylesterase 2; Paraoxonase Nirs; Aromatic esterase 2; Serum arylalkylphosphatase 2; PON2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	354
Samples	Tissue homogenates, cell lysates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.781-50ng/mL
Sensitivity	0.34ng/mL
Intra Assay	CV<10%
Inter Assay	CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	66529294
NCBI GeneID	5445
NCBI Accession #	NP_000296.2
NCBI GenBank Nucleotide #	NM_000305.2
UniProt Accession #	Q15165; O15114; O15115; O75856; Q5FBX7; Q86YL0; A4D1H7; B2RCP9; B4DJD5
Molecular Weight	37,996 Da
NCBI Official Full Name	serum paraoxonase/arylesterase 2 isoform 1
NCBI Official Synonym Full Names	paraoxonase 2

NCBI Official Symbol	PON2
NCBI Protein Information	serum paraoxonase/arylesterase 2
UniProt Protein Name	Serum paraoxonase/arylesterase 2
UniProt Gene Name	PON2
UniProt Synonym Gene Names	PON 2; A-esterase 2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/paraoxonase-2-pon2/66720>